

Dihydralazine sulphate


Catalog No: #S3270

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dihydralazine sulphate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	288.28
Formulation	C ₈ H ₁₂ N ₆ O ₄ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular System
 SMILES: OS(O)(=O)=O.NNC1=NN=C(NN)C2=CC=CC=C12
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting pt: Solubility: Soluble in DMSO
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Stern, H. C. British heart journal. 52 (4): 435C439.

Note: This product is for in vitro research use only