

Metamizole magnesium

Catalog No: #S3272


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Metamizole magnesium
Brief Description	Inhibitors
Purification	99.00%
Target Name	Antibiotic inhibitor
Calculated MW	645
Formulation	C ₂₆ H ₃₂ MgN ₆ O ₈ S ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology
 SMILES: [Mg⁺⁺].CN(CS([O-])(=O)=O)C1=C(C)N(C)N(C1=O)C1=CC=CC=C1.CN(CS([O-])(=O)=O)C1=C(C)N(C)N(C1=O)C1=CC=C1
 Pathways: Microbiology & virology
 Receptor: Antibiotic
 Boiling pt: Melting pt: Solubility: Soluble in Water
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Jage, J. Der Anaesthesist (in German). 57 (4): 382C390.

Note: This product is for in vitro research use only