

Benorilate

Catalog No: #S3276


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Benorilate
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	313.3
Formulation	C17H15NO5
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:ImmunologySMILES:CC(=O)NC1=CC=C(C=C1)OC(=O)C2=CC=CC=C2OC(=O)CPathways:OthersReceptor:OthersBoiling pt:Melting pt:177-181°CSolubility:Soluble in DMSOAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Simila S et al. Eur J Pediatr. 1975 Dec 9;121(1):15-20.

Note: This product is for in vitro research use only