

Tulobuterol hydrochloride

Catalog No: #S3277


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tulobuterol hydrochloride
Brief Description	Inhibitors
Purification	99.00%
Target Name	Adrenergic Receptor agonist
Calculated MW	264.19
Formulation	C ₁₂ H ₁₉ Cl ₂ NO
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: CC(C)(C)NCC(C1=CC=CC=C1Cl)O
 Pathways: GPCR/G Protein Receptor: β 2-adrenergic receptor
 Boiling pt: Melting pt: Solubility: Soluble in DMSO
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.
 Stock solution can be stored below -20 °C for several month.

References

1. Horiguchi, T et al. *Arzneimittelforschung*. 54 (5): 280C285.

Note: This product is for in vitro research use only