

Silver sulfadiazine

Catalog No: #S3279


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Silver sulfadiazine
Brief Description	Inhibitors
Purification	98.00%
Target Name	Antibiotic inhibitor
Calculated MW	357.14
Formulation	C ₁₀ H ₉ AgN ₄ O ₂ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: C1=CN=C(N=C1)[N-]S(=O)(=O)C2=CC=C(C=C2)N.[Ag⁺]
 Pathways: Microbiology & virology Receptor: Antibiotic Binding
 pt: Melting pt: 285°C? Solubility: N/A Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Marx, John et al. Elsevier Health Sciences. p. 814

Note: This product is for in vitro research use only