

Trimeprazine tartrate

Catalog No: #S3282


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Trimeprazine tartrate
Brief Description	Inhibitors
Purification	99.69%
Target Name	Others
Calculated MW	746.98
Formulation	C40H50N4O6S2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: CC(CN1C2=CC=CC=C2SC3=CC=CC=C31)CN(C)C.CC(CN1C2=CC=CC=C2SC3=CC=CC=C31)CN(C)C.[C@@H]([C@H](C(=O)O)O)(C(=O)O)O Pathways: Others Receptor: Others Boiling pt: Melting pt: 153-155°C Solubility: Soluble in DMSO Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Hu OY et al. The Journal of Pharmacy and Pharmacology. 38 (3): 172C6.

Note: This product is for in vitro research use only