

## Cetilistat

Catalog No: #S3283


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Cetilistat
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others inhibitor
Calculated MW	401.58
Formulation	C25H39NO3
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Metabolism systemSMILES:CCCCCCCCCCCCCCCCOC1=NC2=C(C=C(C=C2)C)C(=O)O1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in DMSO, chloroformAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Yamada Y et al. Hormone and Metabolic Research. 40 (8): 539C43.

Note: This product is for in vitro research use only