

Carazolol

Catalog No: #S3284


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Carazolol
Brief Description	Inhibitors
Purification	99.00%
Target Name	Adrenergic Receptor agonist
Calculated MW	298.38
Formulation	C18H22N2O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Cardiovascular SystemSMILES:CC(C)NCC(COC1=CC=CC2=C1C3=CC=CC=C3N2)OPathways:GPCR/G

ProteinReceptor:β-adrenergic receptorBoiling pt:Melting pt:133-137°CSolubility:Soluble in DMSOAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Innis RB et al. Life Sci. 24 (24): 2255C64.

Note: This product is for in vitro research use only