

# Eslicarbazepine Acetate


Catalog No: #S3285

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Eslicarbazepine Acetate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	296.324
Formulation	C17H16N2O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Nervous system  
 SMILES: CC(=O)O[C@H]1CC2=CC=CC=C2N(C3=CC=CC=C3)C(=O)N  
 Pathways: Others  
 Receptor: Others  
 Boiling pt: 183-185°C  
 Melting pt: 183-185°C  
 Solubility: DMSO: 59 mg/mL (199.1 mM)/Ethanol: 11 mg/mL (37.12 mM)  
 Appearance: White Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Sierra-Paredes G, et al. BMC Neurosci. 2014, 15:134.

Note: This product is for in vitro research use only