

Sodium picosulfate monohydrate

Catalog No: #S3286


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sodium picosulfate monohydrate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others inhibitor
Calculated MW	499.42
Formulation	C ₁₈ H ₁₅ NNa ₂ O ₉ S ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: C1=CC=NC(=C1)C(C2=CC=C(C=C2)OS(=O)(=O)[O-])C3=CC=C(C=C3)OS(=O)(=O)[O-].O.[Na+].[Na+] Pathways: Others Receptor: Others
 Boiling pt: Melting pt: Solubility: Soluble in Water Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. FDA approves new colon-cleansing drug for colonoscopy prep

Note: This product is for in vitro research use only