

Octenidine

Catalog No: #S3287


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Octenidine
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	623.82
Formulation	C36H64Cl2N4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: CCCCCCN=C1C=CN(C=C1)CCCCCCCCCN2C=CC(=NCCCCCCC)C=C2.Cl.Cl Pathways: Others Receptor: Others Bol
 ing pt: Melting pt: Solubility: Soluble in DMSO Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Sedlock DM, et al. Microbicidal activity of octenidine hydrochloride, a new alkanediylbis[pyridine] germicidal agent. Antimicrob Agents Chemother. 1985 Dec;28(6):786-90.

Note: This product is for in vitro research use only