


Description

Product Name	Sulfisomidine
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	278.33
Formulation	C12H14N4O2S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:ImmunologySMILES:CC1=CC(=NC(=N1)C)NS(=O)(=O)C2=CC=C(C=C2)NPathways:OthersReceptor:OthersBoiling pt:Melting pt:245CSolubility:Soluble in DMSOAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

N/A

Note: This product is for in vitro research use only