

Amifostine

Catalog No: #S3289


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Amifostine
Brief Description	Inhibitors
Purification	99.20%
Target Name	Others
Calculated MW	214.22
Formulation	C ₅ H ₁₅ N ₂ O ₃ PS
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer SMILES: C(CN)CNCCSP(=O)(O)O Pathways: Others Receptor: Others Boiling pt: Melting pt: 160°C Solubility: Soluble in Water Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Kouvaris JR. Oncologist. 12 (6): 738C47

Note: This product is for in vitro research use only