

## Pikamilone

Catalog No: #S3290


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Pikamilone
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	208.21
Formulation	C10H12N2O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Nervous system  
SMILES: C1=CC(=CN=C1)C(=O)NCCCC(=O)O  
Pathways: Others  
Receptor: Others  
Boiling pt: Melting pt: Solubility: Soluble in DMSO  
Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

N/A

Note: This product is for in vitro research use only