


Description

Product Name	PF-04418948
Brief Description	Inhibitors
Purification	99.00%
Target Name	Prostaglandin Receptor inhibitor
Calculated MW	409.41
Formulation	C23H20FNO5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Endocrine system
SMILES: COC1=CC2=C(C=C1)C=C(C=C2)OCC3(CN(C3)C(=O)C4=CC=C(C=C4)F)C(=O)O
Pathways: GPCR/G Protein Receptor: PGE2
Boiling pt: Melting pt: Solubility: Soluble in DMSO
Appearance: White powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. af Forselles KJ, et al. Br J Pharmacol. 2011 Dec;164(7):1847-1856.

Note: This product is for in vitro research use only