

## Perillyl alcohol

Catalog No: #S3314


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Perillyl alcohol
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	152.23
Formulation	C10H16O
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:CancerSMILES:CC(=C)C1CCC(=CC1)COPathways:OthersReceptor:OthersBoiling pt:228-229°CMelting pt:Solubility:soluble in ethanolAppearance:Viscous oily liquidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Ripple GH. Clin Y Res?2000;6:390-6.

Note: This product is for in vitro research use only