

## Benurestat

Catalog No: #S3315


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Benurestat
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others inhibitor
Calculated MW	228.63
Formulation	C9H9ClN2O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
SMILES: C1=CC(=CC=C1C(=O)NCC(=O)NO)Cl  
Pathways: Others  
Receptor: urease  
Boiling pt: Melting pt: Solubility: Soluble in Ethanol  
Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Andersen JA. Invest Urol. 1975 Mar;12(5):381-6.

Note: This product is for in vitro research use only