

## Fluoromethalone

Catalog No: #S3316


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Fluoromethalone
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	376.46
Formulation	C22H29FO4
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:ImmunologySMILES:C[C@H]1C[C@H]2[C@@H]3CC[C@@]([C@]3(C[C@@H]([C@@]2([C@@]4(C1=CC(=O)C=C4)F)O)C)(C(=O)C)OPath ways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Water: 30mg/L Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Tong Lin, MD. Medicine (Baltimore). 2015 Feb; 94(7): 551.

Note: This product is for in vitro research use only