

Wogonoside

Catalog No: #S3318


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Wogonoside
Brief Description	Inhibitors
Purification	98.66%
Target Name	GABA Receptor inhibitor
Calculated MW	460.39
Formulation	C22H20O11
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system/

Cancer SMILES: COC1=C(C=C(C2=C1OC(=CC2=O)C3=CC=CC=C3)O)O[C@H]4[C@@H]([C@H]([C@@H]([C@H](O4)C(=O)O)O)O)O Pathways: Membrane transporter/Ion channel; Neuroscience Receptor: GABA Binding pt: Melting pt: Solubility: Soluble in DMSO, Methanol Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Hui KM, et al. Biochem. Pharmacol. 2012, 64 (9): 1415C24.

Note: This product is for in vitro research use only