

## Scutellarein

Catalog No: #S3319


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Scutellarein
Brief Description	Inhibitors
Purification	98.31%
Target Name	Src inhibitor
Calculated MW	286.24
Formulation	C15H10O6
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Immunology  
SMILES: C1=CC(=CC=C1C2=CC(=O)C3=C(C(=C(C=C3O)O)O)O)O  
Pathways: Angiogenesis  
Receptor: Src  
Boiling pt: Melting pt: 330°C  
Solubility: DMSO: 38 mg/mL  
Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Sung NY, et al. Korean J Physiol Pharmacol. 2015 Sep;19(5):441-9.

Note: This product is for in vitro research use only