

# (+)-Mk-801 Hydrogen Maleate?


Catalog No: #S3320

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	(+)-Mk-801 Hydrogen Maleate?
Brief Description	Inhibitors
Purification	99.00%
Target Name	NMDAR antagonist
Calculated MW	337.38
Formulation	C <sub>20</sub> H <sub>19</sub> NO <sub>4</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

### Research

Area: Others SMILES: C[C@@]12C3=CC=CC=C3C[C@@H](N1)C4=CC=CC=C24.C(=C\C(=O)O)\C(=O)O Pathways: Neuroscience Receptor: NMDA receptor Boiling pt: Melting pt: Solubility: DMSO: 68 mg/mL (201.55 mM) Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Wong EH, et al. Proc Natl Acad Sci U S A, 1986, 83(18), 7104-7108.

Note: This product is for in vitro research use only