

Oleanonic?Acid

Catalog No: #S3321


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Oleanonic?Acid
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	454.69
Formulation	C30H46O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Immunology/

MetabolismSMILES:C[C@]12CCC(=O)C([C@@H]1CC[C@@]3([C@@H]2CC=C4[C@]3(CC[C@@]5([C@H]4CC(CC5)(C)C)C(=O)O)C)(C)CPathways:OthersReceptor:OthersBoiling pt:551.704CMelting pt:Solubility:Soluble in methanolAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Liu J. Journal of Ethnopharmacology. 1995, 49 (2): 57C68.

Note: This product is for in vitro research use only