

## trans-Chalcone

Catalog No: #S3322


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	trans-Chalcone
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	208.26
Formulation	C15H12O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area:OthersSMILES:C1=CC=C(C=C1)/C=C/C(=O)C2=CC=CC=C2Pathways:OthersReceptor:OthersBoiling pt:208 CMelting pt:55-59  
CSolubility:Soluble in DMSOAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath  
for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Mahapatra. European Journal of Medicinal Chemistry. 92: 839C865.

Note: This product is for in vitro research use only