

Ursonic Acid

Catalog No: #S3323


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ursonic Acid
Brief Description	Inhibitors
Purification	97.00%
Target Name	Others inhibitor
Calculated MW	454.69
Formulation	C30H46O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OthersSMILES:C[C@@H]1CC[C@@]2(CC[C@@]3(C(=CC[C@H]4[C@]3(CC[C@@H]5[C@@]4(CCC(=O)C5(C)C)C)[C@@H]2[C@H]1C)C(=O)O)Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in DMSOAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Yanuar A, et al. Bioinformation. 2014 Feb 19;10(2):52-5.

Note: This product is for in vitro research use only