

## Lycorine

Catalog No: #S3324

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Lycorine
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others inhibitor
Calculated MW	287.31
Formulation	C16H17NO4
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OthersSMILES:C1CN2CC3=CC4=C(C=C3[C@H]5[C@H]2C1=C[C@@H]([C@H]5O)O)OCO4Pathways:OthersReceptor:OthersBoiling pt:Melting pt:253-255°CSolubility:Soluble in DMSOAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Vrijsen R. Biol. Chem. 1986, 261 (2): 505C7.

Note: This product is for in vitro research use only