

## Liquiritigenin

Catalog No: #S3325


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Liquiritigenin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Estrogen/progestogen Receptor agonist
Calculated MW	256.25
Formulation	C15H12O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Others SMILES: C1[C@H](OC2=C(C1=O)C=CC(=C2)O)C3=CC=C(C=C3)O Pathways: Endocrinology/Hormones Receptor: ERβ Binding pt: Melting pt: 206-208°C Solubility: Soluble in DMSO Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Kim, YW. Gastrointestinal and Liver Physiology. 296 (2): G372C81.

Note: This product is for in vitro research use only