

## Angelic Acid


Catalog No: #S3328

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Angelic Acid
Brief Description	Inhibitors
Purification	99.17%
Target Name	Others
Calculated MW	100.12
Formulation	C <sub>5</sub> H <sub>8</sub> O <sub>2</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Immunology  
SMILES: C/C=C(/C)C(=O)O  
Pathways: Others  
Receptor: Others  
Boiling pt: Melting pt: Solubility: Soluble in DMSO  
Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Ernest Small? Culinary herbs, National Research Council Canada, 2006, p. 288

Note: This product is for in vitro research use only