

3,5,4'-Trimethoxystilbene

Catalog No: #S3330


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	3,5,4'-Trimethoxystilbene
Brief Description	Inhibitors
Purification	99.33%
Target Name	Others inhibitor
Calculated MW	270.32
Formulation	C17H18O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:CancerSMILES:COC1=CC=C(C=C1)/C=C/C2=CC(=CC(=C2)OC)OCPathways:OthersReceptor:OthersBoiling pt:Melting pt:57CSolubility:Solubl in ChloroformAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Tsai JH et al. Toxicol Appl Pharmacol. 2013 Nov 1;272(3):746-56

Note: This product is for in vitro research use only