

(1R,2R)-Trans-N-Boc-1,2-Cyclohexanediamine

Catalog No: #S3331


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	(1R,2R)-Trans-N-Boc-1,2-Cyclohexanediamine
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	214.31
Formulation	C11H22N2O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OthersSMILES:CC(C)(C)OC(=O)N[C@@H]1CCCC[C@H]1NPathways:OthersReceptor:OthersBoiling pt:Melting pt:114-115°CSolubility:Soluble in Chloroform, MethanolAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

N/A

Note: This product is for in vitro research use only