

IL-1 β Polyclonal Antibody

Catalog No: #41059

Package Size: #41059-1 50ul #41059-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	IL-1 β Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA IF(paraffin section)
Species Reactivity	Hu Ms Rt
Specificity	IL-1 β Polyclonal Antibody detects endogenous levels of IL-1 β protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human IL-1 β .
Target Name	IL-1 β
Other Names	IL1B; IL1F2; Interleukin-1 beta; IL-1 beta; Catabolin
Accession No.	Swiss-Prot: P01584NCBI Gene ID: 3553
Uniprot	P01584
GeneID	3553;
SDS-PAGE MW	17KD
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/5000.

Not yet tested in other applications.

Images



Western blot analysis of lysate from THP-1/ Mouse brain/Mouse liver , using IL1B Antibody.

Note: This product is for in vitro research use only