

## Cycloastragenol

Catalog No: #S3333


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Cycloastragenol
Brief Description	Inhibitors
Purification	99.62%
Target Name	Telomerase agonist
Calculated MW	490.715
Formulation	C30H50O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Others SMILES: C[C@]12CC[C@@]34C[C@@]35CC[C@@H](C([C@@H]5[C@H])(C[C@H]4[C@@]1(C[C@@H]([C@@H]2[C@@]6(CC[C@@H](O6)C(C)(C)O)C)O)(C)C)OP Pathways: DNA Damage/DNA Repair Receptor: Telomerase Boiling pt: Melting pt: Solubility: Soluble in Chloroform Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

- de Jesus BB, et al. Aging Cell. 2011, 10 (4): 604C21.

Note: This product is for in vitro research use only