

## Sophoridine

Catalog No: #S3339


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Sophoridine
Brief Description	Inhibitors
Purification	99.28%
Target Name	Topoisomerase inhibitor
Calculated MW	248.36
Formulation	C15H24N2O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Cancer SMILES: C1C[C@@H]2[C@H]3CCCN4[C@H]3[C@H](CCC4)CN2C(=O)C1 Pathways: DNA Damage/DNA

Repair Receptor: Topo I Boiling pt: Melting pt: Solubility: Soluble in Chloroform Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Li X, et al. Bioorg Med Chem Lett. 2011, 21(18):5251-4.

Note: This product is for in vitro research use only