

Elemicin


Catalog No: #S3341

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Elemicin
Brief Description	Inhibitors
Purification	97.20%
Target Name	5-HT Receptor agonist
Calculated MW	208.25
Formulation	C12H16O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
SMILES: COC1=CC(=CC(=C1OC)OC)CC=CO
Pathways: GPCR/G Protein; Neuroscience
Receptor: 5-HT2A
Boiling pt: Melting pt: Solubility: Soluble in DMSO
Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Maria Gracia N, et al. Philippine Pharmacists Association, 2014, 3(1)

Note: This product is for in vitro research use only