

Selonsertib?

Catalog No: #S3350


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Selonsertib?
Brief Description	Inhibitors
Purification	95.04%
Target Name	ASK inhibitor
Calculated MW	445.49
Formulation	C24H24FN7O
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:MetabolismSMILES:CC1=C(C=C(C(=C1)F)C(=O)NC2=CC=CC(=N2)C3=NN=CN3C(C)C)N4C=C(N=C4)C5CC5Pathways:ApoptosisReceptor:ASK1Boiling pt:Melting pt:Solubility:DMSO: 36 mg/mL (80.8 mM); Ethanol:24 mg/mL (53.87 mM)Appearance:PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Lin JH, et al. Nephron. 2015, 129(1):29-33.

Note: This product is for in vitro research use only