


Description

Product Name	BIA10-2474
Brief Description	Inhibitors
Purification	99.79%
Target Name	FAAH inhibitor
Calculated MW	300.36
Formulation	C16H20N4O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: O=C(N1C=C(C2=C[N+](O-)=CC=C2)N=C1)N(C3CCCCC3)C Pathways: Metabolism Receptor: FAAH Boiling pt: Melting pt: Solubility: Soluble in DMSO Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. "News Release - Phase I Clinical Trial Rennes". www.bial.com. Retrieved 21 January 2016.

Note: This product is for in vitro research use only