

IL-4R α Polyclonal Antibody

Catalog No: #41061

Package Size: #41061-1 50ul #41061-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IL-4R α Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF ELISA
Species Reactivity	Hu Ms
Specificity	IL-4R α Polyclonal Antibody detects endogenous levels of IL-4R α protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human IL-4R α around the non-phosphorylation site of Y497.
Target Name	IL-4R α
Other Names	IL4R; IL4RA; 582J2.1; Interleukin-4 receptor subunit alpha; IL-4 receptor subunit alpha; IL-4R subunit alpha; IL-4R-alpha; IL-4RA; CD antigen CD124
Accession No.	Swiss-Prot: P24394NCBI Gene ID: 3566
Uniprot	P24394
GeneID	3566;
SDS-PAGE MW	90kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of 293 cells using IL-4Rα Polyclonal Antibody

Note: This product is for in vitro research use only