

## Naloxegol oxalate

Catalog No: #S3355


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Naloxegol oxalate
Brief Description	Inhibitors
Purification	98.00%
Target Name	P450 inhibitor
Calculated MW	741.828
Formulation	C36H55NO15
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Metabolism SMILES: COCCOCCOCCOCCOCCOCCOCCOCCOCCO[C@H]1CC[C@]2([C@H]3CC4=C5[C@]2([C@H]1OC5=C(C=C4)O)CCN3CC=C)O.C(=O)(C(=O)O)OPathways: Metabolism Receptor: CYP3A Boling pt: Melting pt: Solubility: DMSO 40 mg/mL Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Bui K, et al. J Clin Pharmacol. 2016 Aug;56(8):1019-27.

Note: This product is for in vitro research use only