

17 alpha-propionate

Catalog No: #S3357


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-------------------|--|
| Product Name | 17 alpha-propionate |
| Brief Description | Inhibitors |
| Purification | 99.30% |
| Target Name | Androgen Receptor antagonist |
| Calculated MW | 402.52 |
| Formulation | C24H34O5 |
| Storage | 3 years -20°C powder;2 years -80°C in solvent; |

Images



Product Description

Research Area:Endocrine

systemSMILES:CCC(=O)O[C@@]1(CC[C@@H]2[C@@]1(CC[C@H]3[C@H]2CCC4=CC(=O)CC[C@]34C)C(=O)COPathways:Endocrinology/HormonesReceptor:Androgen ReceptorBoiling pt:Melting pt:Solubility:Soluble in DMSOAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Trifu V, et al.. Br J Dermatol. 2011 Jul;165(1):177-83.

Note: This product is for in vitro research use only