

## Cambendazole

Catalog No: #S3361


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Cambendazole
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	302.08
Formulation	C <sub>14</sub> H <sub>14</sub> N <sub>4</sub> O <sub>2</sub> S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Immunology  
 SMILES: CC(C)OC(=O)NC1=CC2=C(C=C1)N=C(N2)C3=CSC=N3  
 Pathways: Others  
 Receptor: Others  
 Boiling pt: Melting pt: Solubility: Soluble in DMSO  
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.  
 Stock solution can be stored below -20 °C for several month.

## References

1. Andrade Sre, et al. 2003 Nov 24;33(4):655-65.

Note: This product is for in vitro research use only