

Eupatilin

Catalog No: #S3362


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Eupatilin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	344.32
Formulation	C18H16O7
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
 SMILES: COC1=C(C=C(C=C1)C2=CC(=O)C3=C(C(=C(C=C3O)O)OC)O)OC
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting pt: Solubility: Soluble in DMSO
 Appearance: Yellow powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.
 Stock solution can be stored below -20 °C for several month.

References

1. Li YY, et al. Oncol Lett. 2015 Oct;10(4):2505-2510.;

Note: This product is for in vitro research use only