

Jervine

Catalog No: #S3363


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Jervine
Brief Description	Inhibitors
Purification	99.00%
Target Name	Hedgehog/Smoothened inhibitor
Calculated MW	425.61
Formulation	C27H39NO3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OthersSMILES:C[C@H]1C[C@@H]2[C@H]([C@H]([C@]3(O2)CC[C@H]4[C@@H]5CC=C6C[C@H](CC[C@@]6([C@H]5C(=O)C4=C3C)C)O)C)NC1Pathways:GPCR/G Protein; Stem CellsReceptor:Hedgehog; SmoothenedBoiling pt:Melting pt:Solubility:Soluble in ChloroformAppearance:White powderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Chen JK, et al. Genes Dev. 2002 Nov 1; 16(21): 2743C2748.;

Note: This product is for in vitro research use only