

Veratramine

Catalog No: #S3364


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Veratramine
Brief Description	Inhibitors
Purification	95.00%
Target Name	Others inhibitor
Calculated MW	409.61
Formulation	C27H39NO2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OthersSMILES:C[C@H]1C[C@H]([C@@H](NC1)[C@@H](C)C2=C(C3=C(C=C2)[C@@H]4CC=C5C[C@H](CC[C@@]5([C@H]4C3)C)O)C)OPa
 thways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in DMSOAppearance:White powderRemark:For obtaining a higher solubility ,
 please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Wang L, et al. Pharmazie. 2008 Aug;63(8):606-10.;2. Li HJ, et al. Chem Pharm Bull (Tokyo). 2006 May;54(5):722-4.;

Note: This product is for in vitro research use only