

## Cephalotaxine

Catalog No: #S3366

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Cephalotaxine
Brief Description	Inhibitors
Purification	95.00%
Target Name	Antiviral inhibitor
Calculated MW	315.37
Formulation	C18H21NO4
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OthersSMILES:COC1=C[C@]23CCCN2CCC4=CC5=C(C=C4[C@@H]3[C@@H]1O)OCO5Pathways:Immunology/InflammationReceptor:AntiviralBoling pt:Melting pt:Solubility:Soluble in Chloroform, DMSOAppearance:White powderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. M A Sells, et al. J Virol. 1988 Aug; 62(8): 2836C2844.;

Note: This product is for in vitro research use only