

Macranthoidin B

Catalog No: #S3367


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-------------------|--|
| Product Name | Macranthoidin B |
| Brief Description | Inhibitors |
| Purification | 99.00% |
| Target Name | Others |
| Calculated MW | 1399.53 |
| Formulation | C65H106O32 |
| Storage | 3 years -20°C powder;2 years -80°C in solvent; |

Images



Product Description

Research

Area:OthersSMILES:C[C@H]1[C@@H]([C@H]([C@H]([C@@H](O1)O[C@@H]2[C@H]([C@H](CO[C@H]2OC3CC[C@@]4[C@H]5CC=C6[C@H]7CC(CC[C@@]7(CC[C@]6([C@@]5(CCC4C3(C)CO)C)C(=O)O[C@H]8[C@@H]([C@H]([C@@H]([C@H](O8)CO[C@H]9[C@@H]([C@H]([C@@H]([C@H](O9)CO)O)O)O)O)O(C)C)O[C@@H]1[C@@H]([C@H]([C@@H]([C@H](O1)CO)O)O)O[C@@H]1[C@H]([C@H]([C@@H]([C@H](O1)CO)O)O)O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in DMSOAppearance:White powderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

- Chen CY, et al.J Chromatogr B Analyt Technol Biomed Life Sci. 2009 Jan 15;877(3):159-65.;

Note: This product is for in vitro research use only