

## Alpha-Hederin

Catalog No: #S3368


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Alpha-Hederin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	750.97
Formulation	C41H66O12
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OthersSMILES:C[C@H]1[C@@H]([C@H]([C@H]([C@@H](O1)O[C@@H]2[C@H]([C@H](CO[C@H]2O[C@H]3CC[C@]4([C@H]([C@]3(C)CO)CC[C@@]5([C@H]4CC=C6[C@]5(CC[C@@]7([C@H]6CC(CC7)(C)C(=O)O)C)C)O)O)O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in Pyridine, MethanolAppearance:White powderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Randhawa MA, et al. Am J Chin Med. 2011;39(6):1075-91.;

Note: This product is for in vitro research use only