

## Hederacoside C

Catalog No: #S3370


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Hederacoside C
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	1221.39
Formulation	C59H96O26
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OthersSMILES:C[C@H]1[C@@H]([C@H]([C@H]([C@@H](O1)O[C@@H]2[C@H](O[C@H]([C@@H]([C@H]2O)O)OC[C@@H]3[C@H]([C@@H]([C@H]([C@@H](O3)OC(=O)[C@@]45CC[C@@]6(C(=CC[C@H]7[C@]6(CC[C@@H]8[C@@]7(CC[C@@H]([C@@]8(C)CO)O[C@H]9[C@@H]([C@H]([C@H](CO9)O)O)O[C@H]1[C@@H]([C@@H]([C@H]([C@@H](O1)C)O)O)C)C)[C@@H]4CC(CC5(C)C)C)O)O)CO)O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:222 CSolubility:Soluble in PyridineAppearance:White powderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Kim JM, et al. Xenobiotica. 2013 Nov;43(11):985-92.;

Note: This product is for in vitro research use only