

Quillaic Acid

Catalog No: #S3371


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Quillaic Acid
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	486.69
Formulation	C30H46O5
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OthersSMILES:C[C@]12CC[C@@H]([C@@]([C@@H]1CC[C@@]3[C@@H]2CC=C4[C@]3(C[C@H]([C@@]5[C@H]4CC(CC5)(C)C)(=O)O)C)(C)C=O)O

Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in DMSOAppearance:White powderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Arrau S, et al. J Ethnopharmacol. 2011 Jan 7;133(1):164-7.;

Note: This product is for in vitro research use only