

Columbianadin

Catalog No: #S3373


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Columbianadin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	328.36
Formulation	C19H20O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: C/C=C(/C)C(=O)OC(C)(C)[C@@H]1CC2=C(O1)C=CC3=C2OC(=O)C=C3 Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Soluble in Chloroform, DMSO Appearance: White powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Zhang YB, et al. Biomedical Chromatography, 2016, 30(2):256-262;

Note: This product is for in vitro research use only