

Dihydrotanshinone I

Catalog No: #S3374


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dihydrotanshinone I
Brief Description	Inhibitors
Purification	95.00%
Target Name	ROS inhibitor
Calculated MW	278.31
Formulation	C18H14O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: CC1COC2=C1C(=O)C(=O)C3=C2C=CC4=C3C=CC=C4 Pathways: Immunology/Inflammation Receptor: ROS
Boiling pt: Melting pt: Solubility: Soluble in Chloroform Appearance: White powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Wang L, et al. *Phytomedicine*. 2015 Nov 15;22(12):1079-87.;

Note: This product is for in vitro research use only