


Description

Product Name	17 α -Estradiol
Brief Description	Inhibitors
Purification	99.00%
Target Name	Estrogen/progestogen Receptor agonist
Calculated MW	272.38
Formulation	C18H24O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Others SMILES: C[C@]12CC[C@H]3[C@H]([C@@H]1CC[C@H]2O)CCC4=C3C=CC(=C4)O Pathways: Endocrinology/Hormones Receptor: ER α ; ER β Boling pt: Melting pt: Solubility: DMSO: 100mM; Ethanol: 100mM Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

N/A

Note: This product is for in vitro research use only